GeneTaq™ 2 × Green Master Mix

Size:

Cat.#	Size
BG10103S	20 R(50 µL vol)
BG10103M	100 R(50 µL vol)
BG10103L	500 R(50 µL vol)

Source: *E.coli* strain with a cloned pol gene from *Thermus aquaticus* YT1.

Storage temperature: -20°C

Exp.: 12 months(-15°C ~ -25°C)

Applications:

1. Routine PCR amplification of DNA fragments up to 5 kb.
2. PCR product for TA cloning.
3. Colony PCR.

www.syngen.tech

Description:

GeneTaq™ 2 × Green Master Mix contains GeneTaq™ DNA Polymerase, dNTPs, other buffer components and stabilizers. GeneTaq™ DNA polymerase is a highly thermostable DNA polymerase, that possesses 5' → 3' polymerase activity, 5' → 3' exonuclease activity, but no detectable 3' → 5' exonuclease (proofreading) activity. GeneTaq™ DNA Polymerase is extremely sensitive and able to amplify trace amounts of DNA in the sample reaction.

GeneTaq™ 2 × Green Master Mix is recommended for any amplification reaction that will be visualized by agarose gel electrophoresis and ethidium bromide staining. It is not recommended when doing any downstream applications in absorbance or fluorescence excitation, as the yellow and blue dyes in the reaction buffer may interfere with these applications.

Reagent supplied:

Cat.#	Size	Quantity	Nuclease-Free Water
BG10103S	20 R(50 µL vol)	0.5 mL	1.5 mL
BG10103M	100 R(50 µL vol)	5 × 0.5 mL	2 × 1.5 mL
BG10103L	500 R(50 µL vol)	25 × 0.5 mL	10 × 1.5 mL

Quality control:

1. Amplification Ability: 25 cycles of PCR amplification of 5 ng Lambda DNA in the expected 5 kb product.
2. Nuclease Activity: No contaminating endonuclease or exonuclease activity has been detected.
3. Protein Purity (SDS-PAGE): GeneTaq™ DNA Polymerase

has >99% pure judged by SDS-PAGE.

4. qPCR for DNA Contamination(*E.coli* Genome): A minimum of 5µL of GeneTaq™ DNA Polymerase is screened for the presence of *E.coli* genomic DNA using qPCR with Taqman Probe specific for *E.coli* 16S rRNA.The measured level of *E.coli* genomic DNA contamination is less than 10 copies of *E.coli* genome.

Protocol:

To prepare several parallel reactions with minimized pipetting errors, prepare a PCR master mix by mixing Nuclease-Free water, primers and GeneTaq™ 2× Master Mix. Prepare sufficient master mix for the number of reactions plus one extra. Aliquot the master mix into individual PCR tubes and then add template DNA.

1. Gently vortex and briefly centrifuge all solutions after thawing.
2. Place a thin-walled PCR tube on ice and add the following component for each 25 µL or 50 µL reaction:

Component	25 µL reaction	50 µL reaction	Final concentration
10 µM forward primer	0.5 µL	1 µL	200 nM
10 µM reverse primer	0.5 µL	1 µL	200 nM
Template DNA	variable	variable	0.1-100 ng
GeneTaq™ 2 × Green Master Mix	12.5 µL	25 µL	1 ×
Nuclease-Free water	to 25 µL	to 50 µL	–

3. Gently vortex the samples and spin down.
4. Perform PCR using recommended thermal cycling

conditions:

Step	Temperature	Time	Number of cycles
Initial Denaturation	94°C	2-5 min	1
Denaturation	94°C	15-30 s	25-40
Annealing	Tm-5°C	10-30 s	
Extension	68°C	10-60 s/kb	1
Final Extension	68°C	5-15 min	

General guidelines:

1. Template:

High quality, purified DNA templates greatly enhances the success of PCR reaction. Recommended amounts of DNA template for a 50 µL reaction are as follows:

DNA	Amount
Genomic DNA	1 ng-1 µg
Plasmid or Viral DNA	1 pg-1 µg

2. Primers:

Oligonucleotide primers are generally 20-40 nucleotides in length and ideally have a GC content of 40-60%. Software such as Oligo 7 can be used to design or analyze primers. The final concentration of each primer in a PCR reaction may be 0.05-1 µM, typically 0.1-0.5 µM.

3. Denaturation:

An initial denaturation of 2 min at 94°C is sufficient for most amplicons from pure DNA templates. For difficult templates such as GC-rich sequences or crude samples, a longer denaturation of 2-8 min at 94°C is recommended

prior to PCR cycling to fully denature the template. With colony PCR, an initial 2-5 min denaturation at 94°C is recommended.

4. Annealing:

The annealing step is typically 15-60 s. Annealing temperature is based on the Tm of the primer pair and is typically 45-68°C. Annealing temperatures can be optimized by doing a temperature gradient PCR starting at 5°C below the calculated Tm.

5. Extension:

The recommended extension temperature is 68°C. Extension time are generally 20-30 s/kb for complex, genomic samples, but can be reduced to 10 s/kb for simple templates (plasmid, virus). To determine specific extension time, we recommend 1 min/kb and reduce the extension time until 10 s/kb. 2 min/kb is recommended for crude samples. A final extension of 5 min at 68°C is recommended.

6. Cycle number:

Generally, 25-35 cycles yields sufficient product. Up to 45 cycles may be required to detect low-copy-number templates.

7. PCR product:

The PCR products generated by GeneTaq™ DNA Polymerase contain dA overhangs at the 3' end, therefore can be ligated to dT/dU-overhang vectors.

References:

1. Chien A, Edgar D B, Trela J M. Deoxyribonucleic-acid polymerase from extreme thermophile *Thermus-aquaticus*[J]. *Journal of Bacteriology*, 1976, 127.

2. Kaledin A S, Sliusarenko A G, Gorodetski S I. Isolation and properties of DNA polymerase from extreme thermophylic bacteria *Thermus aquaticus* YT-1[J]. *Biokhimiia*, 1980, 45(4): 644-651.

Research only, not for therapeutic or diagnostic purposes.
V1.0.2 © 13-14 Beijing SyngenTech Co., Ltd. All rights reserved

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V1.0.2 © 13-14 Beijing SyngenTech Co., Ltd. All rights reserved

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www.syngen.tech

Description:

GeneTaq™ 2 × Colorless Master Mix contains GeneTaq™ DNA Polymerase, dNTPs, other buffer components and stabilizers. GeneTaq™ DNA polymerase is a highly thermostable DNA polymerase, that possesses 5' → 3' polymerase activity, 5' → 3' exonuclease activity, but no detectable 3' → 5' exonuclease (proofreading) activity. GeneTaq™ DNA Polymerase is extremely sensitive and able to amplify trace amounts of DNA in the sample reaction.

Reagent supplied:

Cat.#	Size	Quantity	Nuclease-Free Water
BG10104S	20 R(50 µL vol)	0.5 mL	1.5 mL
BG10104M	100 R(50 µL vol)	5 × 0.5 mL	2 × 1.5 mL
BG10104L	500 R(50 µL vol)	25 × 0.5 mL	10 × 1.5 mL

Quality control:

1. Amplification Ability: 25 cycles of PCR amplification of 5 ng Lambda DNA in the expected 5 kb product.
2. Nuclease Activity: No contaminating endonuclease or exonuclease activity has been detected.
3. Protein Purity (SDS-PAGE): GeneTaq™ DNA Polymerase has >99% pure judged by SDS-PAGE.
4. qPCR for DNA Contamination(*E.coli* Genome): A minimum of 5 µL of GeneTaq™ DNA Polymerase is screened for the presence of *E.coli* genomic DNA using qPCR with Taqman Probe specific for *E.coli* 16S rRNA. The measured level of *E.coli* genomic DNA contamination is less than 10 copies of *E.coli* genome.

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3. Gently vortex the samples and spin down.
4. Perform PCR using recommended thermal cycling conditions:

Step	Temperature	Time	Number of cycles
Initial Denaturation	94°C	2-5 min	1
Denaturation	94°C	15-30 s	25-40
Annealing	Tm-5°C	10-30 s	
Extension	68°C	10-60 s/kb	
Final Extension	68°C	5-15 min	1

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3. Denaturation:

An initial denaturation of 2 min at 94°C is sufficient for most amplicons from pure DNA templates. For difficult templates such as GC-rich sequences or crude samples, a longer denaturation of 2-8 min at 94°C is recommended prior to PCR cycling to fully denature the template. With colony PCR, an initial 2-5 min denaturation at 94°C is recommended.

4. Annealing:

The annealing step is typically 15-60 s. Annealing temperature is based on the Tm of the primer pair and is typically 45-68°C. Annealing temperatures can be optimized by doing a temperature gradient PCR starting at 5°C below the

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Research only, not for therapeutic or diagnostic purposes.

V1.0.2 © 13-14 Beijing SyngenTech Co., Ltd. All rights reserved

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