GeneTaq™ 2 × Colorless Master Mix

Size:

Cat.#	Size
BG10104S	20 R(50 µL vol)
BG10104M	100 R(50 µL vol)
BG10104L	500 R(50 µL vol)

Source: *E.coli* strain with a cloned pol gene from *Thermus aquaticus* YT1.

Storage temperature: -20°C

Exp.: 12 months(-15°C ~ -25°C)

Applications:

- 1. Routine PCR amplification of DNA fragments up to 5 kb.
- 2. PCR product for TA cloning.
- 3. Colony PCR.

www.syngen.tech

Description:

GeneTaq™ 2 × Colorless Master Mix contains GeneTaq™ DNA Polymerase, dNTPs, other buffer components and stabilizers. GeneTaq™ DNA polymerase is a highly thermostable DNA polymerase, that possesses 5' → 3' polymerase activity, 5' → 3' exonuclease activity, but no detectable 3' → 5' exonuclease (proofreading) activity. GeneTaq™ DNA Polymerase is extremely sensitive and able to amplify trace amounts of DNA in the sample reaction.

Reagent supplied:

Cat.#	Size	Quantity	Nuclease-Free Water
BG10104S	20 R(50 µL vol)	0.5 mL	1.5 mL
BG10104M	100 R(50 µL vol)	5 × 0.5 mL	2 × 1.5 mL
BG10104L	500 R(50 µL vol)	25 × 0.5 mL	10 × 1.5 mL

Quality control:

- 1. Amplification Ability: 25 cycles of PCR amplification of 5 ng Lambda DNA in the expected 5 kb product.
- 2. Nuclease Activity: No contaminating endonuclease or exonuclease activity has been detected.
- 3. Protein Purity (SDS-PAGE): GeneTaq™ DNA Polymerase has >99% pure judged by SDS-PAGE.
- 4.qPCR for DNA Contamination(*E.coli* Genome): A minimum of 5µL of GeneTaq™ DNA Polymerase is screened for the presence of *E.coli* genomic DNA using qPCR with Taqman Probe specific for *E.coli* 16S rRNA. The measured level of *E.coli* genomic DNA contamination is less than 10 copies of *E.coli* genome.

Protocol:

To prepare several parallel reactions with minimized pipetting errors, prepare a PCR master mix by mixing Nuclease-Free water, primers and GeneTaq™ 2× Master Mix. Prepare sufficient master mix for the number of reactions plus one extra. Aliquot the master mix into individual PCR tubes and then add template DNA.

- 1. Gently vortex and briefly centrifuge all solutions after thawing.
- 2. Place a thin-walled PCR tube on ice and add the following component for each 25 µL or 50 µL reaction:

Component	25 µL reaction	50 µL reaction	Final concentration
10 µM forward primer	0.5 µL	1 µL	200 nM
10 µM reverse primer	0.5 µL	1 µL	200 nM
Template DNA	variable	variable	0.1-100 ng
GeneTaq™ 2 × Colorless Master Mix	12.5 µL	25 µL	1 ×
Nuclease-Free water	to 25 µL	to 50 µL	–

- 3. Gently vortex the samples and spin down.
- 4.Perform PCR using recommended thermal cycling conditions:

Step	Temperature	Time	Number of cycles
Initial Denaturation	94°C	2-5 min	1
Denaturation	94°C	15-30 s	25-40
Annealing	Tm-5°C	10-30 s	
Extension	68°C	10-60 s/kb	
Final Extension	68°C	5-15 min	1

General guidelines:

- 1. Template:

High quality, purified DNA templates greatly enhances the success of PCR reaction. Recommended amounts of DNA template for a 50 µL reaction are as follows:

DNA	Amount
Genomic DNA	1 ng-1 µg
Plasmid or Viral DNA	1 pg-1 µg

- 2. Primers:

Oligonucleotide primers are generally 20-40 nucleotides in length and ideally have a GC content of 40-60%. Software such as Oligo 7 can be used to design or analyze primers. The final concentration of each primer in a PCR reaction may be 0.05-1 µM, typically 0.1-0.5 µM.

- 3. Denaturation:

An initial denaturation of 2 min at 94°C is sufficient for most amplicons from pure DNA templates. For difficult templates such as GC-rich sequences or crude samples, a longer denaturation of 2-8 min at 94°C is recommended prior to PCR cycling to fully denature the template. With colony PCR, an initial 2-5 min denaturation at 94°C is recommended.

- 4. Annealing:

The annealing step is typically 15-60 s. Annealing temperature is based on the Tm of the primer pair and is typically 45-68°C. Annealing temperatures can be optimized by doing a temperature gradient PCR starting at 5°C below the

calculated Tm.

- 5. Extension:

The recommended extension temperature is 68°C. Extension time are generally 20-30 s/kb for complex, genomic samples, but can be reduced to 10 s/kb for simple templates (plasmid, virus). To determine specific extension time, we recommend 1 min/kb and reduce the extension time until 10 s/kb. 2 min/kb is recommended for crude samples. A final extension of 5 min at 68°C is recommended.

- 6. Cycle number:

Generally, 25-35 cycles yields sufficient product. Up to 45 cycles may be required to detect low-copy-number templates.

- 7. PCR product:

The PCR products generated by GeneTaq™ DNA Polymerase contain dA overhangs at the 3' end, therefore can be ligated to dT/dU-overhang vectors.

References

- 1. Chien A, Edgar D B, Trela J M. Deoxyribonucleic-acid polymerase from extreme thermophile *Thermus-aquaticus*[J]. *Journal of Bacteriology*, 1976, 127.
- 2. Kaledin A S, Sliusarenko A G, Gorodetski S I. Isolation and properties of DNA polymerase from extreme thermophilic bacteria *Thermus aquaticus* YT-1[J]. *Biokhimiia*, 1980, 45(4): 644-651.

Research only, not for therapeutic or diagnostic purposes.

V1.0.2 © 13-14 Beijing SyngenTech Co., Ltd. All rights reserved

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